





Not since we
fitted an . . .



"EXSTAT"
Aerial

No more
interference!

Crystal clear
reception!



The easy and proven way to better radio reception is—install an "EXSTAT" Aerial! In town or country the "EXSTAT" Anti-Interference Aerial will immediately make your radio a better radio.

- VRA 301 Wall Mount Rod only, complete with mounting £4.13.4.
- VRA 302 Chimney Mount Rod only, complete " " £5.13.0.
- ASA 301 Horizontal with Noise Suppression " " £8.12.0.
- ASA 412 Chimney Mount Rod with Noise Suppression £12.17.11.
- ASA 412N Wall Mount Rod with Noise Suppression £12.3.4.
- ASA 304 Extra Aerial Transformers £2.12.0.
- ASA 305 Extra Receiver Transformers £2.12.0.
- ASA 306 Multipoint Transformers £2.15.6.

See your Radio Dealer
about an Exstat Aerial
now!

N.Z. Distributors:
RUSSELL IMPORT CO. LTD.
P.O. Box 102, Wellington.

SEE YOUR SKIN TROUBLE GO IN A FEW DAYS

MAKE THIS 7 DAY TEST

Skin complaints like eczema, spots and rashes don't heal unless you get at their cause — the germs beneath your skin. Valderma, the wonderful new, double-antiseptic balm, gets right under the pores and gives relief almost immediately.

See for yourself how effectively Valderma works! Try this simple 7 day test on your skin.

Rub a little Valderma on your skin where the trouble is.

Do this several times a day for a week.

Within a day or two you will see and feel the difference Valderma makes. Often, in only a few days, your skin will be clear and healthy.



RUB ON



SPOTS GONE

Double Antiseptic Action


Valderma contains two powerful but gentle antiseptics which penetrate deep down under your skin and brings relief to your skin complaint. Its non-greasy, oil-in-water emulsion base, does not clog the pores, and allows septic matter to escape.

Gentle, soothing, healing Valderma does not stain clothing and is invisible on the skin. You simply rub it in, then watch your skin improve. Available everywhere in handy tubes priced at 2/1d.

VALDERMA DOUBLE-ANTISEPTIC
NON-GREASY BALM

Made in New Zealand for Dae Health Laboratories Ltd., London, by Stanford Laboratories Ltd.,
5 Manuka Street, Wellington.

L.I.



RONSONOL

LIGHTER FLUID

Best for all lighters

and

RONSON ACCESSORIES

- Flints 1/- and 2/- per pkt.
- Wired wicks 6d each.
- Literkits 2/- each.
- Service outfits 3/9 each.



2/6

4 oz.
bottle with
special glass
easy-fill pourer.

58.K.10

Ronson Representatives: Ponsford,
Newman & Benson 1948 Ltd.

Ronsonol Distributors:
Delopac (N.Z.) Ltd., Wellington.